

**1 Pilot Marker (triangle drill)**

Mark the bone and carry out initial drilling for directional orientation. Check with paralleling post.

**2 Depth drill 2.5**

Drilling of the final implant length. Use of the drill stop possible. Control with paralleling post.

**3 Final depth drill (based on the implant diameter)**

Final enlargement of the cavity. Drill stops can be used.

**4 Countersink (optional)**

Expansion of the cortical bone.

**5 Thread cutter (manual / automatic)**

Precutting a thread in compact bone (D1/D2).

**6 Manual/machine implant insertion (the SCX-Line as an example)**

Machine insertion: Insert the implant into the bone with the ratchet or the automatic insertion tool.

Manual insertion: Insert the implant into the bone with the ratchet or the automatic insertion tool.

**7 Position locking screw**

Pick up the cover screw with the hexagon screwdriver, insert into the implant and hand-tighten.

Implant diameter (in terms of mechanical capacity)

	3.25	3.75	4.1	4.5	5.5	1	2	3	4	5	6	7
OK	3.25	-	+	-	-	-	-	-	-	-	-	-
	3.75	+	+	+	+	0	-	-	-	-	-	-
	4.1	+	+	+	+	+	+	+	+	+	+	+
	4.5	+	+	+	+	+	+	+	+	+	+	+
UK	5.5	+	+	+	+	+	+	+	+	+	+	+
	5.5	+	+	+	+	+	+	+	+	+	+	+
	4.1	+	+	+	+	+	+	+	+	+	+	+
	3.75	+	+	+	+	0	-	-	-	-	-	-
	3.25	+	+	+	+	-	-	-	-	-	-	-

- = not suitable + = suitable 0 = suitable with reduced masticatory forces below 150 N

Drill sequence with the BEGO Semados® S-Line Tray<sup>Plus</sup>

Implant	Pilot Marker	Depth drill 2.5	Depth drill S 3.25	Depth drill S 3.75	Depth drill S 4.1	Depth drill S 4.5	Depth drill S 5.5	Countersink*	Thread cutter
S-Line 3.25	x	x	x					(x)	(x)
S-Line 3.75	x	x		x				(x)	(x)
S-Line 4.1	x	x		x	x			(x)	(x)
S-Line 4.5	x	x		x	x	x		(x)	(x)
S-Line 5.5	x	x		x	x	x	x	(x)	(x)

x = necessary (x) = optional (depending on bone quality) \*up to marking (see tip of arrow)

Rotational speed

Product	rpm
• Pilot Marker / depth drill	→ 800
• Countersink	→ 800
• Thread cutter	→ 15
• Implant insertion	→ 15

BEGO Implant Systems GmbH & Co. KG  
 Wilhelm-Herbst-Str. 1  
 28359 Bremen, Germany  
 Tel. +49 421 2028-246  
 Fax +49 421 2028-265  
 Email info@bego-implantology.com  
 www.bego.com

The evolution of the BEGO Semados® S-Line

**1993 - 1998**

**1998 - 2017**

**State of the Art**

- Implant shoulder:** Machined (SC implants) or micro-structured (SCX implants) with Platform Switch – depending on the indication and the user preference
- Self-tapping thread:** Increasing core diameter with decreasing thread depth for better primary stability. Progressive microthread for increased bone contact rate
- Optimum cutting grooves:** to better guide the bone graft into the thread base
- Surface:** Trusted high purity and homogeneous TiPure<sup>Plus</sup>-surface

BEGO SEMADOS® S IMPLANTS HAVE BECOME SC AND SCX IMPLANTS.

Your trusted system with an improved design

BEGO Semados® SC/SCX implants: The last word in modern and cost-effective dental implantology

- Quick & functional
- Reliable & proven
- Simple & intuitive
- Universal & successful
- Varied & sophisticated
- 100% German quality product

The BEGO Semados® SC/SCX implant at a glance

- Made of pure grade 4 titanium
- Diameter: 3.25 · 3.75 · 4.1 · 4.5 · 5.5 mm
- Lengths: 7 · 8.5 · 10 · 11.5 · 13 · 15 mm
- Diameter 5.5 in length 7 mm to prevent vertical augmentations
- High purity, homogeneous TiPure<sup>Plus</sup> surface
- Cylindrical implant with rounded apex to protect anatomical structures

www.bego.com/scx

An overview of the SC/SCX implant prosthetics connection

- Platform Switch
- Established conical internal connection with 45° conical taper and hexagon socket for anti-rotation protection
- Inner cone ensures tight fit of the restoration
- CAD/CAM customised prosthetic components made of various materials

REF 84335/02 · CP · © 2016 by BEGO · 2016-11

Partners in Progress

